



ScreenFect[®] siRNA

Transfection Reagent

Product Information & Instruction Manual

Cat. No. of ScreenFect[®] siRNA

S-4001

S-4001-2

ScreenFect



Contents

1. Characteristics of ScreenFect®siRNA	4
2. General Information	4
2.1 Handling of ScreenFect® Transfection Reagents	4
2.2 Storage Conditions	4
2.3 Shipping	5
2.4 Technology	5
3. Basic Considerations for Successful Transfection	5
3.1 Cells	5
3.2 Transfection Methods	5
3.2.1 One-Step Transfection (Combined Plating & Transfection)	5
3.2.2 Two-Step Transfection	6
3.3 Nucleic Acid	6
3.4 ScreenFect® Reagent Dilution and Complex Formation	7
3.5 Dilution Buffer Volumes	7
3.6 Serum and Antibiotics	7
4. ScreenFect®siRNA Transfection Protocols	8
4.1 Overview: How to use ScreenFect®siRNA	9
4.2 siRNA Transfection Protocols	10
4.2.1 Transfection of siRNA in 96-well Plate Format	10
4.2.2 Tabular Protocol for siRNA Transfection in Selected Formats	11
4.3 Optimization of siRNA Transfection	11
4.3.1 Impact of Amounts of siRNA or Reagent	12
on Cell Transfection Results	12
4.3.2 Optimization Protocol for ScreenFect®siRNA Transfection	13
in 96-well Plate Format	13
5. Troubleshooting	16
6. Packaging Sizes with Ordering Information	18



1. Characteristics of ScreenFect®siRNA

ScreenFect®siRNA is a specialized reagent optimized for best performance in siRNA delivery.

ScreenFect®siRNA is serum compatible and free of animal derived components. There is no need for medium change after transfection. ScreenFect® Reagents have relatively low cytotoxicity, allowing easy One-Step¹ cell transfection of recently detached cells. This One-Step procedure reduces the duration of your experiment by one day. Once diluted, ScreenFect®siRNA can be used for a period of up to four days. Always mix the dilution directly before use.

2. General Information

2.1 Handling of ScreenFect® Transfection Reagents

ScreenFect® Reagents should be mixed by vortexing before each use. Do not aliquot and store ScreenFect® Reagents in containers other than the one it is delivered in as contact of undiluted liposomal reagents with plastic surfaces may reduce performance. Dilute by pipetting ScreenFect® Reagents directly into the supplied Dilution Buffer, avoiding contact with the side of tubes, and use pipette action to wash out traces remaining in the pipette tip.

2.2 Storage Conditions

ScreenFect® Reagents should be stored at 2-8°C. Do not leave for extended periods at room temperature and do not freeze. After storage of several months without using the reagent slight precipitation may occur. This has no influence on the performance of ScreenFect® Reagents. If kept refrigerated at 2-8°C the shelflife of ScreenFect® Reagents is guaranteed for at least 12 months.

¹ One-Step transfection: transfection method in which cell plating and transfection is performed in one single step (*see also 3.2.1*)



2.3 Shipping

Shipment is done at ambient temperature to reduce environmental waste and cost. We see no significant loss of activity with storage at room temperature for periods up to two weeks.

2.4 Technology

At ScreenFect, we combine our expertise in chemistry and biology to create a highly interdisciplinary research and development environment that encompasses chemical design, synthesis, liposomal reagent preparation, cell-based screening assays as well as identification and final optimization of novel transfection reagents. With respect to the chemical synthesis, our proprietary thiol-yne based combinatorial click chemistry method allows parallel, high-throughput synthesis of hundreds of novel lipid-like molecules. The method is cost effective, highly efficient and allows for the synthesis of lipids of diverse chemical structures.

3. Basic Considerations for Successful Transfection

3.1 Cells

Cells used for transfection should be in exponential growth phase and have a relatively even density over the entire surface area on which they are plated. We recommend splitting cells when they reach 80 % - 90 % confluence to avoid contact inhibition of cell proliferation. Cells should be mycoplasma free.

3.2 Transfection Methods

Cells can be transfected using the two different methods described below.

3.2.1 One-Step Transfection (Combined Plating & Transfection)

For One-Step transfection (also referred to as **Reverse Cell Transfection**), freshly detached cells in suspension are added to the transfection complexes. The transfection process is thus initiated before cell attachment takes place.



Important facts and benefits of One-Step transfection:

- **Time efficient** procedure (combined plating and transfection)
- **highly recommended for ScreenFect® Reagents**
- Due to the low cytotoxicity of our reagents, reverse cell transfection **does not harm the cells but...**
- **significantly increases transfection efficiencies for most cell lines tested.**




Figure 1: One-Step transfection method

3.2.2 Two-Step Transfection

In the Two-Step or Forward Transfection Method, the cells are plated 24 hours before transfection. The next day, complexes of transfection reagent and nucleic acid are added to the already adherent cells. For more information see *chapter 5. Troubleshooting*.

3.3 Nucleic Acid

For best transfection results, ensure a high quality of your siRNA stock. The amount of siRNA required per transfection is dependent on the particular siRNA construct. Therefore, it is important to determine the appropriate amount of siRNA per transfection through optimization experiments (*see section 4.3*).

Negatively charged nucleic acid (NA) combines with the cationic lipids of ScreenFect® Reagents to form transfection complexes. Transfection complex formation is affected by the ratio of NA to liposomal reagent and has to be optimized (*see section 4.3*).



3.4 ScreenFect® Reagent Dilution and Complex Formation

Transfection complex formation is a critical step for optimal transfection results. The nucleic acid and the transfection reagent must be evenly mixed in Dilution Buffer, both as their separate dilutions as well as when subsequently combined for transfection complex formation. If previously diluted ScreenFect® Reagents are to be used again (up to 4 days storage possible) mix by vortexing immediately before addition to siRNA.

For optimal mixing at the onset of transfection complex formation, equal volumes of nucleic acid dilution and diluted ScreenFect® Reagent are combined using fast pipette action to rapidly form a homogeneous mixture. Vortexing is not recommended during complex formation due to the strong shear forces that may disrupt the complex formation process. For larger scale transfections (e.g. transfection in culture dishes) "splitting" larger volumes into smaller aliquots for the complex formation step is recommended. Ensure at least 20 min of transfection complex formation time for optimal results.

3.5 Dilution Buffer Volumes

To limit unnecessary wastage, we ship 50 mL ScreenFect® Dilution Buffer as standard per mL ScreenFect® reagents. In some situations, however, the amount of ScreenFect® Dilution Buffer may be limiting with respect to the amount of ScreenFect® reagent. If you require additional Dilution Buffer for your particular transfection experiments, please contact us.

3.6 Serum and Antibiotics

Serum does not affect the performance of ScreenFect® Reagents. Although there is no clear evidence for a reduced transfection efficiency using antibiotics, we recommend avoiding Penicillin and Streptomycin during transfection.



4. ScreenFect®siRNA Transfection Protocols

The following pages contain protocols for ScreenFect®siRNA transfection. The reagent volumes suggested in the protocols of *section 4.2* are a guideline. Initial optimization is important. Therefore, an optimization protocol is provided at the end of this chapter. The optimization protocol helps users to determine the optimal amount of ScreenFect®siRNA reagent as well as nucleic acid for their particular transfection experiment and cell line.

We highly recommend the One-Step transfection protocol for all our products. However, for special applications or cell types, the Two-Step protocol may be preferable. Advices on how to use our reagents in Two-Step transfections are given in *section 5 Troubleshooting*.

Note About Dilution Buffer Volumes

To limit unnecessary wastage, we ship 50 mL ScreenFect® Dilution Buffer as standard per mL ScreenFect® reagents, which is sufficient for the vast majority of users. In some situations, however, the amount of ScreenFect® Dilution Buffer may become limiting with respect to the amount of ScreenFect® reagent. If you require additional Dilution Buffer for your particular transfection experiments, please contact us.



4.1 Overview: How to use ScreenFect®siRNA




Figure 2: ScreenFect®siRNA transfection procedure at a glance



4.2 siRNA Transfection Protocols

4.2.1 Transfection of siRNA in 96-well Plate Format

Amounts shown are for one well of a 96-well plate.

1	Dilute 0.25 μl * of ScreenFect®siRNA in Dilution Buffer to a final volume of 7 μl and mix thoroughly. <i>Important: Vortex the reagent once per day of use. Add ScreenFect®siRNA reagent directly into supplied buffer with rapid pipette mixing or vortexing.</i>
2	Dilute 1 pmol* siRNA in the supplied Dilution Buffer to a final volume of 7 μl .
3	Combine the diluted siRNA with the ScreenFect®siRNA dilution and immediately mix using 10 rapid pipette strokes. Allow complex formation to proceed for 20 min at RT. <i>Important: Do not vortex!</i>
4	Add 80 μl freshly detached and resuspended cells to the complexes and mix gently with pipette. <i>Tip: The time-saving one-step cell transfection method may not be suited for all cell types (see section 4.1). To transfect adherent cells, first remove and discard medium from cells, then add 80 μl of fresh culture medium to transfection complexes, mix with pipette and immediately apply to cells.</i>
5	Transfer cells and complexes to one well of a 96-well plate. <i>Note: The time required for efficient endogenous mRNA and protein depletion (24 – 72 h) is gene dependent and has to be determined by the end user.</i>

*Values for amounts of siRNA and reagent given in this table are recommendations. An Optimization Protocol is provided in paragraph 4.3.2. This is intended to help the user to determine the optimal amount of both ScreenFect®siRNA as well as siRNA for transfection of the particular cell type within one experiment.



4.2.2 Tabular Protocol for siRNA Transfection in Selected Formats

For formats other than 96-well, scale reagents using the following table as a guideline.

Format	96-well*	24-well*	6-well*
Step 1	Dilute ScreenFect®siRNA in Dilution Buffer.		
ScreenFect®siRNA	0.25 µl	1 µl	4 µl
Dilution Buffer	7 µl	30 µl	120 µl
Step 2	Dilute siRNA in Dilution Buffer.		
siRNA	1 pmol	5 pmol	25 pmol
Dilution Buffer	7 µl	30 µl	120 µl
Step 3	Add diluted ScreenFect®siRNA to siRNA dilution rapidly and mix with pipette.		
Volume of complexes	14 µl	60 µl	240 µl
Step 4	After 20 min add freshly resuspended cells (or medium) to transfection complexes and transfer mixture to plate.		
Cell suspension (or medium)	80 µl	420 µl	1250 µl

**Values for amounts of siRNA and reagent given in this table are recommendations. An Optimization Protocol is provided in paragraph 4.3.2. This is intended to help the user to determine the optimal amount of both ScreenFect®siRNA as well as siRNA for transfection of the particular cell type within one experiment.*

4.3 Optimization of siRNA Transfection

Optimizing transfection conditions is **highly recommended** for new users to ensure optimal transfection results. Here we provide a basic protocol for varying both the amount of ScreenFect® Reagent as well as siRNA in one convenient experiment.



4.3.1 Impact of Amounts of siRNA or Reagent on Cell Transfection Results

In order to find the best concentrations of siRNA and reagent for your cell line and application, you need to consider the following:

- Higher amounts of siRNA usually result in more efficient gene silencing.
 - Low concentrations of siRNA and reagent result in mild conditions for the cells and therefore in low cell toxicity.
 - Transfection efficiency (→ efficient gene silencing) and cell toxicity negatively correlate with each other.
- It is essential to find working conditions at which these characteristics are balanced. Be sure to choose robust conditions ("optimal working range" in the picture below) that are not sensitive to experimental deviations.




Figure 3: General correlations between reagent volumes, transfection efficiency and cell viability.

Note: Efficiency and toxicity effects are dependent on the cell type used.



4.3.2 Optimization Protocol for ScreenFect®siRNA Transfection in 96-well Plate Format

1.	Prepare 3 different siRNA dilutions:																
1A	Take a total of 4 pmol, 8 pmol or 12 pmol siRNA and dilute each amount in Dilution Buffer to a final volume of 40 µl.																
1B	Take 3 x 10 µl from <i>each</i> of the above three siRNA samples and distribute between the wells of a 96-well PCR tube plate as shown in the pipetting scheme under 3. (10 µl of diluted reagent will be added to these nine samples later.)																
	<i>Note: The siRNA samples can alternatively be distributed between tubes of PCR-strips.</i>																
2.	Prepare three different dilutions of ScreenFect®siRNA transfection reagent:																
	We recommend to use PCR stripes for dilution preparation.																
	<table><tr><th>Sample ID</th><th>Corresponding final volume of ScreenFect® Reagent per 96-well:</th><th>Volume of Dilution Buffer:</th><th>Volume of ScreenFect®siRNA or :</th></tr><tr><td>2.1</td><td>0.15 µl</td><td>197 µl</td><td>3 µl</td></tr><tr><td>2.2</td><td>0.25 µl</td><td>78 µl</td><td>2 µl</td></tr><tr><td>2.3</td><td>0.35 µl</td><td>193 µl</td><td>7 µl</td></tr></table>	Sample ID	Corresponding final volume of ScreenFect® Reagent per 96-well:	Volume of Dilution Buffer:	Volume of ScreenFect®siRNA or :	2.1	0.15 µl	197 µl	3 µl	2.2	0.25 µl	78 µl	2 µl	2.3	0.35 µl	193 µl	7 µl
Sample ID	Corresponding final volume of ScreenFect® Reagent per 96-well:	Volume of Dilution Buffer:	Volume of ScreenFect®siRNA or :														
2.1	0.15 µl	197 µl	3 µl														
2.2	0.25 µl	78 µl	2 µl														
2.3	0.35 µl	193 µl	7 µl														
	<i>Important: Add ScreenFect® Reagents directly into supplied buffer and mix with rapid pipette action or vortexing.</i>																

Note: slightly higher volumes of Dilution Buffer are used in the optimization protocol to allow greater accuracy in pipetting small volumes, which can be beneficial for less experienced users.



3.	Mix each diluted ScreenFect® sample with the different siRNA samples:																																																																																																																																																		
	Add 10 µl of each diluted ScreenFect®siRNA sample to each of the three different amounts of siRNA in a PCR-reaction tube, mix and incubate 20 min at RT.																																																																																																																																																		
	Take 10 µl aliquots from the three ScreenFect® Reagent dilutions and add to the siRNA samples containing 1 pmol (samples 1-3 for plate coordinates C5-C7 as shown below). <i>Immediately mix with 10 rapid pipette strokes, do not vortex!</i> Repeat for the siRNA samples containing 2 pmol (samples 4-6) and 3 pmol (samples 7-9).																																																																																																																																																		
	Pipetting scheme:																																																																																																																																																		
	<table><tr><th></th><th>1</th><th>2</th><th>3</th><th>4</th><th>5</th><th>6</th><th>7</th><th>8</th><th>9</th><th>10</th><th>11</th><th>12</th><th rowspan="9">siRNA amount: 1 pmol 2 pmol 3 pmol</th></tr><tr><th>A</th><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr><tr><th>B</th><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr><tr><th>C</th><td></td><td></td><td></td><td></td><td>1</td><td>2</td><td>3</td><td></td><td></td><td></td><td></td><td></td></tr><tr><th>D</th><td></td><td></td><td></td><td></td><td>4</td><td>5</td><td>6</td><td></td><td></td><td></td><td></td><td></td></tr><tr><th>E</th><td></td><td></td><td></td><td></td><td>7</td><td>8</td><td>9</td><td></td><td></td><td></td><td></td><td></td></tr><tr><th>F</th><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr><tr><th>G</th><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr><tr><th>H</th><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr><tr><td colspan="5">µl of ScreenFect®siRNA</td><td>0.15</td><td>0.25</td><td>0.35</td><td colspan="5"></td><td></td></tr><tr><td colspan="5">per well (and sample ID):</td><td>2.1</td><td>2.2</td><td>2.3</td><td colspan="5"></td><td></td></tr></table>		1	2	3	4	5	6	7	8	9	10	11	12	siRNA amount: 1 pmol 2 pmol 3 pmol	A													B													C					1	2	3						D					4	5	6						E					7	8	9						F													G													H													µl of ScreenFect®siRNA					0.15	0.25	0.35							per well (and sample ID):					2.1	2.2	2.3						
	1	2	3	4	5	6	7	8	9	10	11	12	siRNA amount: 1 pmol 2 pmol 3 pmol																																																																																																																																						
A																																																																																																																																																			
B																																																																																																																																																			
C					1	2	3																																																																																																																																												
D					4	5	6																																																																																																																																												
E					7	8	9																																																																																																																																												
F																																																																																																																																																			
G																																																																																																																																																			
H																																																																																																																																																			
µl of ScreenFect®siRNA					0.15	0.25	0.35																																																																																																																																												
per well (and sample ID):					2.1	2.2	2.3																																																																																																																																												
4.	Mix transfection complexes with cells:																																																																																																																																																		
	The One-Step transfection method gives best results. Add 80 µl freshly dissociated cells to the complexes and then mix gently with pipette before transferring to a 96-well cell culture plate.																																																																																																																																																		
	<i>For the Two-Step method, first remove the used medium from the 96-well cell culture plate of already adherent cells. Then add 80µl medium to transfection complexes and transfer medium and complexes to plate with adherent cells.</i>																																																																																																																																																		

Note: The optimization protocol (Steps 1-4) can be downloaded as a single paged PDF from our homepage. Information regarding the optimization of siRNA transfections in different formats is provided on the next page.



4.3.2.1 Optimization of siRNA Transfections in Selected Formats

The following table is intended to assist the user with the optimization of siRNA transfections in different formats. Only the volume ranges of reagent, Dilution Buffer and the ranges of nucleic acid amounts are included into the table below. For operational instructions, please refer to *section 4.3.2*.

Format	96-well	24-well	6-well
Step 1	Dilute ScreenFect® Reagent in Dilution Buffer.		
ScreenFect® Reagent	0.15, 0.25, 0.35 µl	0.75, 1, 1.5 µl	3, 4, 5 µl
Total volume of dilution	10 µl	30 µl	120 µl
Step 2	Dilute nucleic acid in Dilution Buffer.		
Nucleic acid	1, 2 pmol	5, 10 pmol	25, 40 pmol
Total volume of dilution	10 µl	30 µl	120 µl
Step 3	Add diluted ScreenFect® Reagent rapidly to the nucleic acid dilutions and mix with pipette.		
Volume of complexes	20 µl	60 µl	240 µl
Step 4	After 20 min add freshly resuspended cells (or medium in case of Two-Step method) to complexes and transfer mixture to cell culture plate.		
Cell suspension	80 µl	420 µl	1250 µl

Note: slightly higher volumes of Dilution Buffer are used for some of the optimization protocols. The reason for this is to allow greater accuracy in pipetting small volumes.



5. Troubleshooting

In case of low transfection efficiency, use the troubleshooting guide below as a basis to identify the problem.

What might be the problem in case of low transfection efficiency?

Possible causes

Suggested solutions

Poor quality of siRNA or insufficient siRNA amount

- Ensure a high quality and the right concentration of your siRNA stock solution.
- Optimise the concentration of siRNA according to the optimization protocol in *section 4.3.2*.

Insufficient complex formation

- Make sure the diluted ScreenFect® Reagent and siRNA solutions are, after combination, mixed immediately using rapid pipette action. Do not vortex! (*See section 3.4*)
- If the volume used for complex formation is greater than about 200 µl, it may help to split the complex formation step into several tubes of reduced volumes – in order to make the initial mixing step more efficient (*see section 3.4*).
- Incubate the prepared mix of siRNA and ScreenFect® Reagent for at least 20 minutes at RT.

Improper ratio of siRNA to reagent

- Optimize the ratio of siRNA to ScreenFect®siRNA according to the optimization protocol in *section 4.3.2*.

Incorrect storage

- ScreenFect® Reagents should be stored at 2–8°C. Do not leave for extended periods at room temperature and do not freeze. Try to avoid excessive warming-cooling if used frequently by users. Remove from 2–8°C for brief periods when needed and replace at 2–8°C.

Cells are unhealthy

- The cells used for transfection should be in exponential growth phase and have even density over surface area. Make sure your cells are free of any contamination. The use of antibiotics is recommended during passaging. Ensure the density of your cell culture does not get too high or too low during the experiment and while passaging cells – always maintain timely passaging.
-



Inappropriate cell density
or cells are not actively
dividing

- The cells used for transfection should be in exponential growth phase and have even density over the entire surface area. Ensure appropriate density of your cell culture by timely passaging. This is important for both One- and Two-Step transfection protocols.
- For the Two-Step method, split your cell culture early enough before transfection (15–24 h in advance).
- Excessive passaging of cells decreases the transfection performance so use cells with similar passage number between experiments to ensure reproducibility.
- Use a new batch of cells.

Inhibition of complex
formation

- Serum does not affect the performance of ScreenFect® Reagents, but we recommend avoiding antibiotics during transfection as well as anionic inhibitors such as EDTA or dextran sulphate.

Improper protocol

- We highly recommend the One-Step transfection method for all our reagents. Based on our experience, One-Step transfection leads to increased performance of ScreenFect® Reagents compared to Two-Step.
 - The time-saving One-Step transfection method may not be suited for all cell types and applications. Only if One-Step transfection does not lead to the desired results test the Two-Step method: Cells are plated into the desired format the day before so that they are adherent at time of transfection. Shortly before the addition of transfection complexes to adherent cells, we recommend removing all (for weakly adherent cells, half) of the medium from the cells. In case of 96-well plates, 80 µl of fresh medium is then mixed with the transfection complexes and this mixture is then added to the adherent cells. Take care to avoid cell detachment, especially for weakly adherent cells.
-



6. Packaging Sizes with Ordering Information

ScreenFect®siRNA

is provided in the following packaging sizes:

- 0.2 ml ScreenFect®siRNA + 10 ml Dilution Buffer (Cat. No.: S-3001-2)
- 1 ml ScreenFect®siRNA + 50 ml Dilution Buffer (Cat. No.: S-3001)

For questions about ScreenFect® Reagents or optimization of transfection protocols or in case of any problems with our products, please feel free to contact us:

info@screenfect.com

Or visit our website for more information:

www.screenfect.com



ScreenFect GmbH
KIT Innovation Center
Hermann-von-Helmholtz-Platz 6
76344 Eggenstein-Leopoldshafen
Germany

Fax +49 (0)7247 3983012

info@screenfect.com

www.screenfect.com